

User Manual for

m r M L M

multi-locus random-SNP-effect Mixed Linear Model tools for
multi-locus GWAS and linkage analysis

(version 2.0)

**Ren Wen-Long, Ni Yuan-Li, Wen Yang-Jun, Wang Shi-Bo,
Huang Bo, Zhang Yuan-Ming (soyzzhang@mail.hzau.edu.cn)**

Last updated on December 18, 2016

Disclaimer: While extensive testing has been performed by Yuan-Ming Zhang's Lab (Statistical Genomics Lab) at Huazhong Agricultural University and Nanjing Agricultural University, the results are, in general, reliable, correct or appropriate. However, results are not guaranteed for any specific datasets. We strongly recommend that users validate the mrMLM results with other software packages, such as EMMAX, GAPIT v2 and PLINK.


Download website:

<https://cran.r-project.org/web/packages/mrMLM/index.html>

Citation:

Method	Paper
mrMLM	Wang et al. <i>Scientific Reports</i> 2016, 6 :19444 ^[1]
GCIM	Wang et al. <i>Scientific Reports</i> 2016, 6 :29951 ^[2]
FASTmrEMMA	Wen et al. <i>Briefings in Bioinformatics</i> 2016, Accepted ^[3]

Note: These references are listed in section of Reference.



This work was supported by the National Natural Science Foundation of China (31571268), Huazhong Agricultural University Scientific & Technological Self-innovation Foundation (Program No. 2014RC020) and Nanjing Agricultural University.



Contents

1 INTRODUCTION	4
1.1 Why mrMLM ?	4
1.2 Getting started	4
1.2.1 Online install	5
1.2.2 Offline install	5
1.2.2.1 Install GTK+	5
1.2.2.2 Install the add-on packages	5
1.2.2.3 Install mrMLM	6
2 DATASETS	6
2.1 GWAS datasets	7
2.1.1 Phenotypic data	7
2.1.2 Genotypic data	7
2.1.2.1 Numeric format	7
2.1.2.2 Character format	8
2.1.2.3 Hapmap format	9
2.1.3 Kinship	9
2.1.4 Population structure	10
2.2 Linkage analysis datasets	11
2.2.1 GCIM format	11
2.2.1.1 Phenotypic data	11
2.2.1.2 Genotypic data	12
2.2.1.3 Linkage map data	12
2.2.2 WinQTLCart format	13
2.2.3 QTL IciMapping format	14
2.2.3.1 Phenotypic data	14
2.2.3.2 Genotypic data	14
2.2.3.3 Linkage map data	15

2.2.4 Covariate format	15
3 ANALYSIS AND RESULTS	16
3.1 mrMLM module	16
3.1.1 Input data	16
3.1.2 Run program	18
3.1.3 Output results	20
3.1.4 Manhattan plot	22
3.1.5 QQ plot	23
3.2 FASTmrEMMA module	25
3.2.1 Input data	25
3.2.2 Run program	28
3.2.3 Output results	29
3.2.4 Manhattan plot	31
3.2.5 QQ plot	32
3.3 GCIM module	33
3.3.1 Input data and Parameter setting	34
3.3.2 Run program	36
3.3.3 Output results	38
3.3.4 Draw plot	39
4 REFERENCE	41

1 INTRODUCTION

1.1 Why mrMLM?

mrMLM (**m**ulti-locus **r**andom-SNP-effect **M**ixed **L**inear **M**odel) software package is an R package with interactive graphic user interface (GUI) for multi-locus genome-wide association study (GWAS) and multi-locus linkage analysis. At present it (version 2.0) includes three modules, which are 1) mrMLM, 2) FASTmrMLM (**F**AST **m**ulti-locus **r**andom-SNP-effect **E**MMMA), and 3) GCIM (**G**enome-wide **C**omposite **I**nterval **M**apping). The first two modules (mrMLM and FASTmrMLM) are used to conduct multi-locus GWAS, the last one (GCIM) is used to implement multi-locus linkage analysis.

(i) **mrMLM**, aims to provide a user-friendly interface to conduct multi-locus GWAS via mrMLM methodology and to visualize its results. Its visualization is based on package **qqman**, which can help draw figures such as Manhattan and QQ plots.

(ii) **FASTmrEMMA**, aims to provide a user-friendly interface to conduct multi-locus GWAS via FASTmrEMMA methodology and to visualize its results. Its visualization is also based on package **qqman**, which can help draw figures such as Manhattan and QQ plots.

(iii) **GCIM**, aims to provide a user-friendly interface to conduct multi-locus linkage analysis via GCIM methodology and to visualize its results. Its visualization is based on package **graphics**, which can help draw figures about the plot of LOD score (y_1 axis) and $-\log_{10}(\text{P-value})$ (y_2 axis) against genome position (cM).

mrMLM 2.0 is able to work on the popular platforms, like Windows, Linux (desktop) and MacOS. And the interactive GUI is based on available add-on package **RGtk2**, via the aid of another package **gWidgetsRGtk2**.

1.2 Getting started

mrMLM is a package that runs in the R software environment, which can be freely downloaded from <https://cran.r-project.org/web/packages/mrMLM/index.html>, or request from the maintainer, Dr Yuan-Ming Zhang (soyzhang@mail.hzau.edu.cn or soyzhang@hotmail.com).

1.2.1 Online install

Within R environment, the mrMLM software can be installed directly using the below command:

```
install.packages(pkgs="mrMLM")
```

1.2.2 Offline install

1.2.2.1 Install GTK+

You may need to install GTK+ before installing RGtk2, because RGtk2 depends on GTK+.

For **Windows** user, you do as below:

Download GTK+ here

(<http://sourceforge.net/projects/gladewin32/files/gtk%2B-win32-runtime/2.10.11/gtk-2.10.11-win32-1.exe>).

Run the resulting file ([gtk-2.10.11-win32-1.exe](#)), which is an automated installer that will help you to complete the installation of Gtk2 libraries. If you use **64-bit** R software, please download corresponding GTK+ version.

For **Mac OS** users, you do as below:

Download GTK+ here (<http://sourceforge.net/projects/gtk-osx/files/latest/download>).

Extract and run the resulting file ([gtk-osx-docbook-1.2.tar.gz](#)).

For **Linux** users, you do as below:

You may or may not upgrade the GTK libraries depending on your distribution.

There are more details on RGtk2 at RGtk2's home page (<http://www.ggobi.org/rgtk2/>).

1.2.2.2 Install the add-on packages

The following R packages are needed: RGtk2, cairoDevice, gWidgets, gWidgetsRGtk2, RGtk2Extras and qqman, which can be downloaded from [CRAN](https://cran.r-project.org/) (<https://cran.r-project.org/>). To install them in order, as some depend on others. Within R environment, these packages can be installed directly using the below command:

```
install.packages(pkgs=c("RGtk2","cairoDevice","gWidgets","gWidgetsRGtk2","RGtk2Extras", "qqman","openxlsx","stringr"))
```

1.2.2.3 Install mrMLM

Open R GUI, select "Packages"—"Install package(s) from local files..." and then find the mrMLM package which you have downloaded on your desktop.

Within R environment, launch the mrMLM by command: `library(mrMLM)`, then the following dialog will appear.




Figure 1.1. Screenshot of mrMLM package GUI

To restart the GUI, the command `mrMLM()` can be issued.

Note: Users should install RStudio software (<https://www.rstudio.com/>) in order to obtain the [User Manual.pdf](#) file from the **User Manual** button. If the RStudio package is not installed, users can decompress the mrMLM package and find the User Manual file (name: [Instruction.pdf](#)) in the folder of ".../mrMLM/inst/doc".

2 DATASETS

2.1 GWAS datasets

2.1.1 Phenotypic data

The **Phenotypic** file should be a ***.csv** format file. The first column stands for individual ID, such as 33-16, Nov-38 and 4226. The second column is phenotypic values for the trait. Note that the phenotypic file includes only one trait. The first element in the first column must be "<Phenotype>" or "<Trait>".

	A	B
1	<Phenotype>	
2	B46	42
3	B52	72.5
4	B57	41
5	B64	74.5
6	B68	65
7	B73	83.25
8	B73HTRHM	73
9	B75	56
10	B76	53
11	B77	50
12	B79	66.5
13	CM7	28
14	CML10	110
15	CML103	88.5
16	CML108	55.5

Figure 2.1. The Phenotypic file

2.1.2 Genotypic data

The **Genotypic** file should be a ***.csv** format file.

2.1.2.1 Numeric format

The first column, named "**rs#**" in the first row, stands for marker ID. The second column, named "**chrom**" in the first row, stands for chromosome. The third column, named "**pos**" in the first row, stands for the position (bp) of SNP in the chromosome. The fourth column, named "**genotype for code 1**" in the first row, stands for reference bases. If the base is missing, the observed base for this marker in the first individual is what we list. And each of the remaining columns is for one individual. In their first rows, the individual names are appeared. For each marker, homozygous genotypes are expressed by 1 and -1, respectively, and the heterozygous and missing genotypes are indicated by zero. Note that the genotypes with code **1** will be listed in the **Result** files. For the [FASTmrEMMA](#) module, homozygous genotypes are expressed by 1 (AA) and 0 (aa), respectively, and the heterozygous and missing genotypes are indicated by 0.5.

	A	B	C	D	E	F	G	H
1	rs#	chrom	pos	genotype	33-16	Nov-38	A4226	A4722
2	PZB00859.1	1	157104	C	1	1	1	1
3	PZA01271.1	1	1947984	C	1	-1	1	-1
4	PZA03613.2	1	2914066	G	1	1	1	1
5	PZA03613.1	1	2914171	T	1	1	1	1
6	PZA03614.2	1	2915078	G	1	1	1	1
7	PZA03614.1	1	2915242	T	1	1	1	1
8	PZA02117.1	1	223466480	A	1	1	1	-1
9	PZA00403.5	1	223466873	T	1	1	1	0
10	PZB01979.2	1	224421551	A	1	-1	1	-1
11	PZB01979.3	1	224421962	C	1	1	-1	1
12	PZA00658.21	1	224518409	T	0	1	1	-1

Figure 2.2 The genotypic file with numeric format in mrMLM module

	A	B	C	D	E	F	G	H
1	rs#	chrom	pos	genotype	33-16	Nov-38	A4226	A4722
2	PZB00859.1	1	157104	C	1	1	1	1
3	PZA01271.1	1	1947984	C	1	0	1	0
4	PZA03613.2	1	2914066	G	1	1	1	1
5	PZA03613.1	1	2914171	T	1	1	1	1
6	PZA03614.2	1	2915078	G	1	1	1	1
7	PZA03614.1	1	2915242	T	1	1	1	1
8	PZA02117.1	1	223466480	A	1	1	1	0
9	PZA00403.5	1	223466873	T	1	1	1	0.5
10	PZB01979.2	1	224421551	A	1	0	1	0
11	PZB01979.3	1	224421962	C	1	1	0	1
12	PZA00658.21	1	224518409	T	0.5	1	1	0

Figure 2.3 The genotypic file with numeric format in FASTmrEMMA module

2.1.2.2 Character format

The first three columns are same as those in the “**Figure 2.2 The genotypic file with numeric format in mrMLM module**”. The differences are that the marker values are character, such as **A**, **T**, **C**, **G** and **N**, and the other notations are heterozygous genotypes. The “**N**” indicates missing. The first rows from the fourth to last columns are individual code.

	A	B	C	D	E	F	G	H
1	rs#	alleles	chrom	pos	33-16	Nov-38	A4226	A4722
2	PZB00859.1	A/C	1	157104	C	C	C	C
3	PZA01271.1	C/G	1	1947984	C	G	C	G
4	PZA03613.2	G/T	1	2914066	G	G	G	G
5	PZA03613.1	A/T	1	2914171	T	T	T	T
6	PZA03614.2	A/G	1	2915078	G	G	G	G
7	PZA03614.1	A/T	1	2915242	T	T	T	T
8	PZA02117.1	A/G	1	2.23E+08	A	A	A	G
9	PZA00403.5	C/T	1	2.23E+08	T	T	T	N
10	PZB01979.2	A/G	1	2.24E+08	A	G	A	G
11	PZB01979.3	C/G	1	2.24E+08	C	C	G	C
12	PZA00658.21	C/T	1	2.25E+08	N	T	T	C

Figure 2.4 The genotypic file with character format

2.1.2.3 Hapmap format

Please see the TASSEL software in details. Here we introduce simply. The first eleven columns describe the specific information of markers and individuals, and their column names must be "rs#", "alleles", "chrom", "pos", "strand", "assembly#", "center", "protLSID", "assayLSID", "panel" and "QCcode". In the "rs#" (1), "chrom" (3) and "pos" (4) columns, their information is described as the above. The values for marker genotypes should be character, such as **AA**, **TT**, **CC**, **GG**, **NN**, **AC** and **AG**, where the "NN" indicates missing or unknown genotypes. In the 2 and 5 to 11 columns, the **no available** information must be marked by "NA". All the individual genotypic information will be showed from the 12 to last columns, the first element in each column is individual ID (name) and the others are the genotypes (character).

	A	B	C	D	E	F	G	H	I	J	K	L
1	rs#	alleles	chrom	pos	strand	assembly#	center	protLSID	assayLSID	panel	QCcode	33-16
2	PZB00859.1	A/C	1	157104	+	AGPv1	Panzea	NA	NA	maize282	NA	CC
3	PZA01271.1	C/G	1	1947984	+	AGPv1	Panzea	NA	NA	maize282	NA	CC
4	PZA03613.2	G/T	1	2914066	+	AGPv1	Panzea	NA	NA	maize282	NA	GG
5	PZA03613.1	A/T	1	2914171	+	AGPv1	Panzea	NA	NA	maize282	NA	TT
6	PZA03614.2	A/G	1	2915078	+	AGPv1	Panzea	NA	NA	maize282	NA	GG
7	PZA03614.1	A/T	1	2915242	+	AGPv1	Panzea	NA	NA	maize282	NA	TT
8	PZA02117.1	A/G	1	223466480	+	AGPv1	Panzea	NA	NA	maize282	NA	AA
9	PZA00403.5	C/T	1	223466873	+	AGPv1	Panzea	NA	NA	maize282	NA	TT
10	PZB01979.2	A/G	1	224421551	+	AGPv1	Panzea	NA	NA	maize282	NA	AA
11	PZB01979.3	C/G	1	224421962	+	AGPv1	Panzea	NA	NA	maize282	NA	CC
12	PZA00658.21	C/T	1	224518409	+	AGPv1	Panzea	NA	NA	maize282	NA	NN

Figure 2.5 The genotypic file with Hapmap format

2.1.3 Kinship

The **Kinship** file should be a *.csv file. The number of rows (or columns) equals to the number of the common individuals between the phenotypic and genotypic datasets. If the Kinship matrix is calculated by this software, we calculate only the Kinship

matrix between the common individuals. If the Kinship matrix has been obtained and uploaded from a known file, it is possible that the number and order of individuals in the known Kinship file are not consistent with those of the common (valid) individuals in the Phenotype and Genotype files. At this situation, the software will change the known K matrix in order that the number and order of new K matrix matches the number and order of common (valid) individuals in the Phenotypic and Genotypic files. In the known K matrix, the **first element** in the first column must be **the number of valid individuals**, such as 263; and the other elements in the first column are individual ID (names). If the number of markers is very large, i.e., 50,000, we recommend that users calculate the K matrix using the other softwares, especially for FASTmrEMMA.

	A	B	C	D	E	F
1	263					
2	33-16	1.00809	0.45954	0.50677	0.42503	0.45591
3	Nov-38	0.45954	1.03352	0.43048	0.47044	0.39597
4	A4226	0.50677	0.43048	1.01717	0.45409	0.43775
5	A4722	0.42503	0.47044	0.45409	0.89002	0.34874
6	A188	0.45591	0.39597	0.43775	0.34874	1.0099
7	A214N	0.34693	0.33421	0.39779	0.29244	0.33058
8	A239	0.43593	0.46499	0.40323	0.36691	0.39597
9	A272	0.34874	0.40505	0.31423	0.3887	0.44138
10	A441-5	0.47952	0.44138	0.47226	0.47952	0.49224
11	A554	0.39779	0.45954	0.5431	0.48679	0.4214
12	A556	0.50858	0.40505	0.45954	0.40142	0.40687

Figure 2.6. The Kinship file

2.1.4 Population structure

The **Population Structure** file is a *.csv file. Using the **Structure** software, the population structure matrix may be calculated. The first column stands for the valid individual ID (names). The first and second elements in the first column must be "<Covariate>" and "<Trait>" respectively. If population structure matrix has k columns, please input all the k columns. In the second row, it must be "Q1", "Q2", ... , "Q k " following the "<Trait>".

If the **Population Structure** file has been obtained and uploaded from a known file, it is possible that the number and order of individuals in the known file are not consistent with those of the common (valid) individuals in the further analysis. At this

situation, the software will change the known matrix in order that the number and order of new **Population Structure** matrix match the number and order of common (valid) individuals in the Phenotypic and Genotypic files.

	A	B	C
1	<Covariate>		
2	<Trait>	Q1	Q2
3	33-16	0.972	0.014
4	Nov-38	0.993	0.004
5	A4226	0.917	0.012
6	A4722	0.854	0.111
7	A188	0.982	0.005
8	A214N	0.017	0.221
9	A239	0.963	0.002
10	A272	0.122	0.859
11	A441-5	0.531	0.464
12	A554	0.979	0.002

Figure 2.7. The Population Structure file

2.2 Linkage analysis datasets

2.2.1 GCIM format

2.2.1.1 Phenotypic data

The **Phenotypic** file should be a *.csv format file. In the first column, “phenotype” is appeared in the first row and **individual ID or individual names** are listed in the other rows, such as DH6-10, DH6-101 and DH6-101. In the second column, the phenotypic values for the first trait are listed. The trait name is appeared in the first row, such as DS1-BLUEs, and observations are listed in the other rows. In the other columns, it is similar to the column 2. Missing phenotypic values is indicated by **NA**.

	A	B
1	phenotype	T19
2	DH6-10	75.33
3	DH6-101	105
4	DH6-102	96.33
5	DH6-104	81
6	DH6-105	101.67
7	DH6-108	89
8	DH6-11	82.67
9	DH6-111	133.33
10	DH6-112	105
11	DH6-114	NA
12	DH6-119	104.67
13	DH6-124	75.67
14	DH6-125	113.33

Figure 2.8. The Phenotypic file with GCIM format

2.2.1.2 Genotypic data

The **Genotypic** file should be a *.csv format file. In the first column, "**genotype**" is showed in the first row and marker names are appeared from row 2 to row $m+1$ (m is the number of markers). In the column 2, **individual ID or individual names** is appeared in the first row and marker genotypes for all the markers are listed in the other rows. For the other rows, it is similar to the column 2.

Population Type. At present, GCIM can analyze five population types derived from two parental lines. $F_1 = P_1 \times P_2$.

1. BC₁: $F_1 \times P_1$. 1 stands for AA, -1 stands for Aa, and 99 stands for missing genotypes.
2. BC₂: $F_1 \times P_2$. 1 stands for Aa, -1 stands for aa, and 99 stands for missing genotypes.
3. DH: doubled haploid lines derived from F_1 . 1 stands for AA, -1 stands for aa, and 99 stands for missing genotypes.
4. RIL: recombinant inbred lines derived from repeatedly selfing since F_1 . 1 stands for AA, -1 stands for aa, and 99 stands for missing genotypes.
5. Chromosome Segment Substitution Line (CSSL): 1 stands for AA, -1 stands for aa, 99 stands for missing genotypes.

	A	B	C	D	E	F	G	H
1	genotype	DH6-10	DH6-101	DH6-102	DH6-104	DH6-105	DH6-108	DH6-11
2	RGA3(1)	-1	99	-1	1	-1	-1	1
3	wPt-6358	-1	99	99	99	99	-1	1
4	Hplc2	1	1	-1	1	-1	-1	1
5	wPt-9752	1	99	99	99	99	99	1
6	abc156a	1	1	-1	1	-1	-1	-1
7	RGA36b(2)	1	99	-1	1	99	-1	-1
8	bcd98	-1	1	-1	1	-1	-1	-1
9	wmc24	-1	1	-1	1	-1	-1	-1
10	ksuG9c	-1	1	-1	1	-1	-1	-1
11	wPt-2436	-1	99	99	99	99	-1	-1
12	wPt-4886	-1	99	99	99	99	-1	-1

Figure 2.9. The genotypic file with GCIM format

2.2.1.3 Linkage map data

The **Linkage Map** file should be a *.csv format file. In the first column, “marker” is appeared in the first row, and marker names are listed in the other rows. In the second column, “chr” is chromosome and appeared in the first row, and chromosome information for all the markers are listed in the other rows. In the third column, “pos” is marker position information, marker position (cM) are listed in the other rows. Note that the marker position is not interval length.

	A	B	C
1	marker	chr	pos
2	RGA3(1)	1	0
3	wPt-6358	1	3.034
4	Hplc2	1	8.8291
5	wPt-9752	1	10.1452
6	abc156a	1	41.3408
7	RGA36b(2)	1	43.8429
8	bcd98	1	51.9122
9	wmc24	1	54.5814
10	ksuG9c	1	55.3333
11	wPt-2436	1	59.1871
12	wPt-4886	1	61.7209
13	wmc120	1	62.4012
14	cdo105	1	63.0679

Figure 2.10. The linkage map file with GCIM format

2.2.2 WinQTLCart format

Please see the WinQTLCart software in details. Here we introduce simply. Between the first “–start” and “–stop” information, marker information is listed, including **chromosome name, marker name and position (interval length)**. Between the second “–start” and “–stop”, genotypic information is showed. Between the third “–start” and “–stop”, phenotypic information is appeared. If variable information is included, it will be **otrait** information between the fourth “–start” and “–stop”.

```
#FileID 1438801284
#bychromosome
-type position
-function 1
-Units cM
-chromosomes 3
-maximum 13
-named yes
-start
-Chromosome C1
AXR-1      0.0000
HH.335C-COL 9.3000
EC.480C    17.2000
EC.650C    29.9000
GC.845L    38.7000
G23554     52.8000
CH.160L-COL 57.8000
DF.701L    72.4000
AD.333C    76.6000
GBB.80L    93.2000
GD.97L     97.0000
FD.8411L   115.5000
MXF-F-001  116.5000
```

Figure 2.11. The mcd file in the WinQTLCart software

2.2.3 QTL IciMapping format

2.2.3.1 Phenotypic data

All the information for each trait is listed in one row. The first column stands for trait ID. And each of the remaining columns stands for the phenotypic value of every individual. Note that -100 stands for missing phenotypes.

	A	B	C	D	E	F	G	H	I	J	K	L
1	T15	103	111.7	106.7	92.67	122	93.67	108	118.3	102.3	-100	130
2	T19	75.33	105	96.33	81	101.7	89	82.67	133.3	105	-100	104.7

Figure 2.12. The phenotypic format in the QTL IciMapping software

2.2.3.2 Genotypic data

All the information for each marker is listed in one row. The first column stands for marker ID. And each of the remaining columns stands for one individual. In QTL IciMapping, 2 was used to represent the marker genotype of the first parent (P_1), 0 for the second parent (P_2), 1 for the F_1 , and -1 for missing genotypes.

	A	B	C	D	E	F	G	H	I	J	K
1	RGA3(1)	0	-1	0	2	0	0	2	2	2	-1
2	wPt-6358	0	-1	-1	-1	-1	0	2	2	2	-1
3	Hplc2	2	2	0	2	0	0	2	2	2	0
4	wPt-9752	2	-1	-1	-1	-1	-1	2	2	2	-1
5	abc156a	2	2	0	2	0	0	0	0	2	0
6	RGA36b(2)	2	-1	0	2	-1	0	0	0	2	-1
7	bcd98	0	2	0	2	0	0	0	0	2	0
8	wmc24	0	2	0	2	0	0	0	0	2	0
9	ksuG9c	0	2	0	2	0	0	0	0	2	0
10	wPt-2436	0	-1	-1	-1	-1	0	0	0	0	-1
11	wPt-4886	0	-1	-1	-1	-1	0	0	0	0	-1
12	wmc120	0	2	0	2	0	0	0	0	0	2

Figure 2.13. The genotypic file in the QTL IciMapping software

2.2.3.3 Linkage map data

The linkage map file has three columns. The first column stands for marker ID, the second column stands for chromosome, and the third column is marker position (cM) in the chromosome, which is not the length of marker interval.

	A	B	C
1	RGA3(1)	1	0
2	wPt-6358	1	3.034
3	Hplc2	1	8.8291
4	wPt-9752	1	10.1452
5	abc156a	1	41.3408
6	RGA36b(2)	1	43.8429
7	bcd98	1	51.9122
8	wmc24	1	54.5814
9	ksuG9c	1	55.3333
10	wPt-2436	1	59.1871
11	wPt-4886	1	61.7209
12	wmc120	1	62.4012

Figure 2.14. The Linkage map format in the QTL IciMapping software

2.2.4 Covariate format

The **Covariate** file should be a *.csv format file. In the first column, “covariate” is appeared in the first row and **individual names** or **individual ID**, such as DH6-10, DH6-101 and DH6-101, are listed in the other rows. In the second column, the name for the first covariate is appeared in the first row and covariate information is listed in the other row. The other columns are similar to the second column if there are several covariates.

	A	B
1	covariate	EnvironName
2	DH6-10	A
3	DH6-101	A
4	DH6-102	A
5	DH6-104	A
6	DH6-105	A
7	DH6-108	A
8	DH6-11	A
9	DH6-111	A
10	DH6-112	A
11	DH6-114	A
12	DH6-119	A

Figure 2.15. The covariate file with GCIM format

3 ANALYSIS AND RESULTS

3.1 mrMLM module

Click the button “[Multi-locus Random-SNP-effect Mixed Linear Model \(mrMLM\)](#)”, then the following dialog will appear.




Figure 3.1.1. Screenshot of mrMLM module GUI

3.1.1 Input dataset

Use the **Input Dataset** button to input dataset files, and then a dialog box will be appeared. In the dialog box, there are four steps. First, users select the dataset formats, which include **mrMLM numeric** format, **mrMLM character** format and **hapmap**

format used in the TASSEL software. Then, use the **Genotype** and **Phenotype** buttons to input the genotypic and phenotypic datasets, respectively. Once one file is successfully uploaded, one tabbed page is added to the software notebook. Third, two things will be implemented in this step. One is to sort the individuals between the genotypic and phenotypic files and all the common individuals between the two files are selected to be analyzed in the further analyses. Another is to transfer the character genotypes into the numeric genotypes if the genotypes are character. Once users press the **DO** button, the two things will be conducted. Once the two files will be successfully uploaded, two tabbed pages (Genotype and Phenotype) will be added to the software notebook. Finally, use the **Kinship** and **Population Structure** buttons to input the kinship and population structure matrices, respectively. If one file is successfully uploaded, the corresponding data page will be added to the notebook. Note that the **Kinship** and **Population Structure** buttons have two options. For the **kinship** button, one is to directly upload the kinship matrix and another is to calculate the kinship matrix using this software. For the **Population Structure** button, the population structure matrix may be not included in the mixed linear model of the GWAS if it has no effect on GWAS. If not, it should be included in the mixed model. The population structure matrix in your uploaded file will be deleted one column if the sum of all the Q-matrix columns for one individual (one row) equal to 1. In the filter dialog, you should choose one column that should be deleted.




Figure 3.1.2. The Input Dataset dialog of mrMLM module




Figure 3.1.3. The Kinship dialog of mrMLM module




Figure 3.1.4. The population structure dialog of mrMLM module




Figure 3.1.5. Filter dialog of the population structure of mrMLM module

3.1.2 Run program

Use the **Parameter Setting** button to set parameters before running the program. “Search radius of candidate gene (kb)” means to keep the one marker with having the least P-value, and to delete all the other markers within the radius of the associated marker with the least P-value. Use the **Run** button to execute the software. If the program runs, a progress bar with the “**Please be patient...**” words will appear in the bottom of the interface. If the program finished, a bar with the “**All done.**” will appear.

Parameter Setting

1. Critical P-value in rMLM:

0.01

2. Search radius of candidate gene (kb):

20

3. Critical LOD score in mrMLM:

3

OK Cancel

Figure 3.1.6. The Parameter Setting dialog of mrMLM module

Multi-locus Random-SNP-effect Mixed Linear Model (mrMLM)

Input Datasets

Parameter Settings

Run

Save

Manhattan Plot

QQ Plot

User Manual

Exit

About the software

	V1	V2	V3	V4	V5	V6	V7	V8	V9	V10	V11	V12	V13	V14	V15
1	263														
2	33-16	1.00809	0.45954	0.50677	0.42503	0.45591	0.34693	0.43593	0.34874	0.47952	0.39779	0.50858	0.43775	0.43048	0.3887
3	Nov-38	0.45954	1.03352	0.43048	0.47044	0.39597	0.33421	0.46499	0.40505	0.44138	0.45954	0.40505	0.42503	0.36509	0.46499
4	4226	0.50677	0.43048	1.01717	0.45409	0.43775	0.39779	0.40323	0.31423	0.47226	0.5431	0.45954	0.45228	0.4323	0.46136
5	4722	0.42503	0.47044	0.45409	0.89002	0.34874	0.29244	0.36691	0.3887	0.47952	0.48679	0.40142	0.38507	0.37781	0.32695
6	A188	0.45591	0.39597	0.43775	0.34874	1.0099	0.33058	0.39597	0.44138	0.49224	0.4214	0.40687	0.37781	0.36509	0.39415
7	A214N	0.34693	0.33421	0.39779	0.29244	0.33058	1.0208	0.36509	0.37054	0.34693	0.39052	0.36327	0.32513	0.29607	0.68114
8	A239	0.43593	0.46499	0.40323	0.36691	0.39597	0.36509	1.02443	0.40142	0.45591	0.43775	0.45409	0.45046	0.31242	0.41777
9	A272	0.34874	0.40505	0.31423	0.3887	0.44138	0.37054	0.40142	0.96086	0.49042	0.45954	0.37962	0.43048	0.34693	0.39779
10	A441-5	0.47952	0.44138	0.47226	0.47952	0.49224	0.34693	0.45591	0.49042	1.02988	0.45409	0.45409	0.44501	0.34693	0.39415
11	A554	0.39779	0.45954	0.5431	0.48679	0.4214	0.39052	0.43775	0.45954	0.45409	1.01717	0.40868	0.40505	0.41958	0.46862
12	A556	0.50858	0.40505	0.45954	0.40142	0.40687	0.36327	0.45409	0.37962	0.45409	0.40868	0.99901	0.46136	0.37599	0.39234
13	A6	0.43775	0.42503	0.45228	0.38507	0.37781	0.32513	0.45046	0.43048	0.44501	0.40505	0.46136	1.0208	0.32513	0.41777
14	A619	0.43048	0.36509	0.4323	0.37781	0.36509	0.29607	0.31242	0.34693	0.34693	0.41958	0.37599	0.32513	1.01717	0.34148
15	A632	0.3887	0.46499	0.46136	0.32695	0.39415	0.68114	0.41777	0.39779	0.39415	0.46862	0.39234	0.41777	0.34148	1.01535
16	A634	0.36509	0.42322	0.43956	0.29244	0.41232	0.73926	0.37599	0.35783	0.36146	0.46318	0.36509	0.36691	0.31968	0.8991
17	A635	0.40323	0.42866	0.45409	0.3106	0.41413	0.67024	0.43048	0.35056	0.39415	0.41232	0.37417	0.37236	0.33603	0.81737
18	A654	0.40323	0.47044	0.49224	0.39234	0.38144	0.33603	0.4995	0.37962	0.45409	0.57579	0.42503	0.45046	0.31787	0.46318
19	A659	0.48497	0.42685	0.40868	0.37599	0.3996	0.41413	0.45591	0.37417	0.36327	0.41777	0.4105	0.35419	0.51585	0.43956
20	A661	0.36327	0.36691	0.33603	0.38325	0.34511	0.30515	0.36327	0.38144	0.33785	0.4105	0.37417	0.31787	0.37236	0.35056
21	A679	0.29607	0.49769	0.30697	0.33421	0.30152	0.39597	0.26337	0.36327	0.33785	0.43593	0.34329	0.33785	0.36691	0.47952
22	A680	0.25793	0.45046	0.30515	0.27791	0.27427	0.46862	0.25429	0.3324	0.31605	0.38144	0.31242	0.32695	0.31605	0.52675
23	A682	0.41232	0.37054	0.44864	0.45046	0.41595	0.34511	0.35964	0.31242	0.39234	0.3996	0.41232	0.30333	0.38507	0.35238

Please be patient ...

Figure 3.1.7. A running program interface of mrMLM module




Figure 3.1.8. A finished program interface (the rMLM Results: [Result1](#))




Figure 3.1.9. A finished program interface (the mrMLM Results: [Result2](#))

3.1.3 Output results

Use **Save** button to save the results as *.csv files. The **Results in the rMLM** are saved as [Result1.csv](#) and the **Results in the mrMLM** are saved as [Result2.csv](#). If click **OK** button (after the **Save** button), a dialog is used to choose the pathway and the file name for the saving files.




Figure 3.1.10. The result **Save** dialog for the rMLM and mrMLM methods




Figure 3.1.11. The Save dialog

The **Result1** table with twelve columns shows the results from the rMLM (random-SNP-effect mixed linear model) method. The corresponding column names are as follows: reference sequence number (rs#, marker name), chromosome, marker's position (bp) in the chromosome, population mean value (Mean), residual variance (σ^2 , Sigma2), priori variance of the k th SNP effect (ϕ_k^2 , Sigma2_k), SNP effect (γ_k , Effect), posteriori variance of SNP effect ($\text{var}(\gamma_k)$, Sigma2_k_posteriori), Wald test statistic value, the P-value of Wald test, significance and genotype for code 1, respectively. In the significance column, only significant markers under the critical value $0.05/m_e$ are marked.

	A	B	C	D	E	F	G	H	I	J	K	L
1	RS#	Chromosome	Marker Position (bp)	Mean	Sigma2	Sigma2_k	SNP effect	Sigma2_k_posteriori	Wald	P_vald	Genotype for code 1	Significance
2	PZB00859.1	1	157104	59.455	179.6497	7.6588	2.3981	1.3757	3.0389	0.0813	C	
3	PZA01271.1	1	1947984	60.8254	180.1372	4.8869	-1.8452	1.2174	2.2973	0.1296	C	
4	PZA03613.2	1	2914066	60.9935	181.7164	0.0082	0.0011	0.0907	1.00E-04	0.9905	C	
5	PZA03613.1	1	2914171	60.993	181.7165	0.0082	-8.34E-06	0.0907	8.45E-09	0.9999	T	
6	PZA03614.2	1	2915078	60.9925	181.7154	0.0082	0.0035	0.0907	0.0015	0.9688	C	
7	PZA03614.1	1	2915242	60.9927	181.7161	0.0082	0.0023	0.0907	6.00E-04	0.9801	T	
8	PZA02117.1	1	223466480	60.9935	181.7164	0.0082	0.0013	0.0907	2.00E-04	0.9888	A	
9	PZA00403.5	1	223466873	60.9933	181.7162	0.0082	0.0021	0.0907	5.00E-04	0.9819	T	
10	PZB01979.2	1	224421551	60.9912	181.7143	0.0082	-0.0051	0.0907	0.0032	0.955	A	
11	PZB01979.3	1	224421962	60.9921	181.7149	0.0082	0.0044	0.0907	0.0023	0.9615	C	
12	PZA00658.21	1	224518409	60.993	181.7165	0.0082	6.93E-05	0.0907	5.84E-07	0.9994	T	
13	PZA00658.19	1	224518576	60.9918	181.716	0.0082	0.0026	0.0907	8.00E-04	0.9775	A	
14	FRM2478.22	1	226781494	61.3264	181.3825	2.0867	-0.8241	1.185	0.4836	0.4868	C	

Figure 3.1.12. Results in the rMLM (Result1)

The **Result2** table with seven columns shows the final results of the mrMLM (multi-locus random-SNP-effect mixed linear model) method. The corresponding column names are as follows: reference sequence number (rs#, marker names), chromosome, marker's position (bp) in the chromosome, QTN effect, LOD score, the proportion of phenotypic variance explained by the putative QTN, and genotype for code 1, respectively.

	A	B	C	D	E	F	G	H	I
1	RS#	Chromosome	Marker Position (bp)	QTN effect	LOD score	r ² (%)	Genotype	for code 1	
2	PZA03188.4	1	280719882	8.3065	6.5787	10.5893	C		
3	PZA03559.1	2	15810363	5.4929	4.2506	6.7945	C		
4	PZA00112.5	5	13664679	-8.1972	4.8377	7.1934	A		
5									
6									
7									
8									

Figure 3.1.13. Results in the mrMLM (Result2)

3.1.4 Manhattan plot

Use **Manhattan plot** button to preview Manhattan plot in an independent dialog window. Before saving the Figure, please set the width and height of the Figure, with the unit of pixel (px). And set word resolution in the Figure, with the unit of 1/72 inch, being pixels per inch (ppi). And set figure resolution in the Figure, with the unit of pixels per inch (ppi). The colors of the two adjacent chromosomes can be changed via the combo box, with a drop-down option. Set the critical value for $-\log_{10}(\text{P-value})$, which is defaulted the value of $0.05/m_e$, where m_e is the effective number of markers (please see Wang et al. *Scientific Reports* 2016, 6: 19444). Use **Save** button to choose a path and to save the Figure, with three frequently used image formats: *.png, *.tiff and *.jpeg.




Figure 3.1.14. The Manhattan Plot of mrMLM module using example data

3.1.5 QQ plot

Use **QQ plot** button to preview the QQ plot in an independent dialog window. The settings of the width, height, word resolution and figure resolution are the same as those in the Manhattan plot. Set the critical P-value for QQ plot, which is defaulted the value of 0.95. This is because the P-values are a mixture of a χ^2 distribution with one degree of freedom and a point mass at one. Note that users may also change this 0.95 based on yourself results.




Figure 3.1.15. The QQ Plot of mrMLM module using example data

Using [example data in this software](#), the preview Manhattan and QQ plots with the general resolution parameter settings were shown in **Figures 3.1.14** and **3.1.15**, respectively. Using [Arabidopsis real data in Atwell et al. \(2010\) Nature 465: 627-631](#), the preview Manhattan and QQ plots with the general resolution parameter settings were shown in **Figures 3.1.16** and **3.1.17**, respectively.




Figure 3.1.16. The Manhattan Plot of mrMLM module using Arabidopsis real data




Figure 3.1.17. The Q-Q Plot of mrMLM module using Arabidopsis real data

If users want to obtain **high resolution** figure, we recommend the parameter settings shown in **Figure 3.1.18** and **Figure 3.1.19** for Manhattan plot and QQ plot, respectively. It is emphasized that preview figure may be not the same with the **saved figure**.

Width (px):	<input type="text" value="10000"/>	Height (px):	<input type="text" value="6000"/>
Word resolution (1/72 inch, ppi):	<input type="text" value="30"/>	Figure resolution (ppi):	<input type="text" value="300"/>
Chromosome color (odd):	<input type="text" value="blue"/>	Chromosome color (even):	<input type="text" value="red"/>
Critical value for Manhattan Plot:	<input type="text" value="3.2314"/>	<input type="button" value="Save"/> <input type="button" value="Cancel"/>	

Note: The parameter settings of preview figure are for general resolution, see manual for parameter settings for high resolution figure (*.png,*.tiff,*.jpeg).

Figure 3.1.18. Parameter settings for high resolution Manhattan plot of mrMLM Module

Width (px):	<input type="text" value="10000"/>	Height (px):	<input type="text" value="6000"/>
Word resolution (1/72 inch, ppi):	<input type="text" value="30"/>	Figure resolution (ppi):	<input type="text" value="300"/>
Point color:	<input type="text" value="black"/>	Line color:	<input type="text" value="red"/>
Critical P-value of deleting points:	<input type="text" value="0.95"/>	<input type="button" value="Save"/> <input type="button" value="Cancel"/>	

Note: The parameter settings of preview figure is for general resolution, see manual for parameter settings of high resolution figure (*.png,*.tiff,*.jpeg).

Figure 3.1.19. Parameter settings for high resolution QQ plot of mrMLM module

3.2 FASTmrEMMA module

Click the button “[Fast Multi-locus Random-SNP-effect EMMA \(FASTmrEMMA\)](#)”, then the following dialog will appear.




Figure 3.2.1. Screenshot of FASTmrEMMA module GUI

3.2.1 Input Data

Use the **Input Dataset** button to input dataset files, and then a dialog box will be appeared. In the dialog box, there are five steps. First, users select the dataset formats, which include **FASTmrEMMA numeric** format, **FASTmrEMMA character** format

and **hapmap** format used in the TASSEL software. Second, use the **Genotype** and **Phenotype** buttons to input the genotypic and phenotypic datasets, respectively. Once one file is successfully uploaded, one tabbed page is added to the software notebook. Third, users select one objective function: **Restricted Likelihood Function or Likelihood Function**. Fourth, some things will be implemented in this step: 1) firstly, to delete missing phenotypic data, then to sort the individuals between the genotypic and phenotypic files and all the common individuals between the two files are selected to be analyzed in the further analyses; 2) to transfer the character genotypes into the numeric genotypes if the genotypes are character; and 3) to select one objective function. Once users press the **DO** button, these things will be implemented. Once the two files will be successfully uploaded, two tabbed pages (Genotype and Phenotype) will be added to the software notebook. Finally, use the **Kinship** and **Population Structure** buttons to input the kinship and population structure matrices, respectively. If one file is successfully uploaded, the corresponding data page will be added to the notebook. Note that there are two options for the **Kinship** and **Population Structure** buttons. For the **kinship** button, one is to directly upload the kinship matrix and another is to calculate the kinship matrix in this software. For the **Population Structure** button, the population structure matrix is not included in FASTmrEMMA methodology if it has no effect on GWAS. If not, it should be included in FASTmrEMMA methodology. The population structure matrix in your uploaded file will be deleted one column if the sum of all the Q-matrix columns for one individual (one row) equal to 1. In the filter dialog, you should choose one column that should be deleted.

Note: If the number of SNP markers in genotypic dataset is large (number of SNP markers is more than 100,000), please input the kinship matrix file directly. About the **input file format**, please see **Direction 1** in the end of the manual.




Figure 3.2.2. The Input Dataset dialog of FASTmrEMMA module




Figure 3.2.3. The Kinship dialog of FASTmrEMMA module




Figure 3.2.4. The population structure dialog of FASTmrEMMA module




Figure 3.2.5. Filter dialog of the population structure of FASTmrEMMA module

3.2.2 Run Program

Use the **Parameter Setting** button to set parameters before run the program. “Search radius of candidate gene (kb)” means to keep the one marker with having the least P-value, and to delete all the other markers within the radius of the associated marker with the least P-value. Use the **Run** button to execute the software. If the program runs, a progress bar with the “**Please be patient...**” words will appear in the bottom of the interface. If the program finished, a bar with the “**All done.**” will appear.




Figure 3.2.6. The Parameter Setting dialog of FASTmrEMMA module




Figure 3.2.7. A running program interface of FASTmrEMMA module




Figure 3.2.8. A finished program interface of FASTmrEMMA module (**Results: Result1**)




Figure 3.2.9. A finished program interface of FASTmrEMMA module (**Results: Result1**)

3.2.3 Output results

Use **Save** button to save the results as *.csv files. The **Result1** are saved as **Result1.csv** and the **Result2** are saved as **Result2.csv**. If click **OK** button (after the **Save** button), a dialog is used to choose the pathway and the file name for the saving

files.

Note: About the **explanation of Result1 and Result2** in details, please see **Direction 2** in the end of the manual.




Figure 3.2.10. The result **Save** dialog of FASTmrEMMA module




Figure 3.2.11. The **Save** dialog of FASTmrEMMA module

The **Result1** table with twelve columns shows the results. The corresponding column names are as follows: **Marker** (marker name or reference sequence number), **Chromosome**, **Marker position (bp)** in the chromosome, the **P_value** for each QTN, **SNP effect**, **QTN-to-residual variance ratio** (ratio of QTN variance to residual variance), **QTN variance**, **Residual variance** for the current SNP, **MAF** (minor allele frequency), **r² (%)** (the proportion of phenotypic variance explained by the putative

QTN), **genotype for code 1**, and **Critical P-value** (critical value for the significant test, $0.05/m_e$), respectively.

	A	B	C	D	E	F	G	H	I	J	K	L	M
1	Marker	Chromosome	Marker position (bp)	F-value	SNP effect	QTN-to-residual variance ratio	QTN variance	Residual variance	MAF	r ² (%)	Genotype for code 1	Critical F-value	
2	PZB00859.1	1	157104	0.1161	4.553	0.1625	28.2184	173.6541	0.2015	0.8062	C	1.00E-04	
3	PZA01271.1	1	1947984	0.1654	-3.4311	0.1001	17.4347	174.1106	0.4639	0.6752	C		
4	PZA03613.2	1	2914066	1	8.00E-04	4.54E-05	0.008	175.5027	0.289	3.18E-08	G		
5	PZA03613.1	1	2914171	1	-0.00017597	4.54E-05	0.008	175.5028	0.2319	1.32E-09	T		
6	PZA03614.2	1	2915078	1	0.0021	4.54E-05	0.008	175.5024	0.4601	2.59E-07	G		
7	PZA03614.1	1	2915242	1	0.0015	4.54E-05	0.008	175.5026	0.4335	1.18E-07	T		
8	PZA02117.1	1	223466480	1	0.001	4.54E-05	0.008	175.5027	0.3802	5.26E-08	A		
9	PZA00403.5	1	223466873	1	0.0013	4.54E-05	0.008	175.5026	0.3118	8.96E-08	T		
10	PZB01979.3	1	224421551	1	-0.0024055	4.54E-05	0.008	175.5023	0.346	3.07E-07	A		
11	PZB01979.3	1	224421962	1	0.002	4.54E-05	0.008	175.5024	0.2548	1.88E-07	C		
12	PZA00658.21	1	224518409	1	2.00E-04	4.54E-05	0.008	175.5028	0.1483	1.69E-09	T		
13	PZA00658.19	1	224518576	1	0.0012	4.54E-05	0.008	175.5027	0.1977	5.28E-08	A		
14	PHM2478.22	1	226781494	0.3718	-1.8025	0.053	9.2878	175.1404	0.1369	0.0948	C		
15	PZA00664.3	1	227542648	0.3973	-1.0066	0.0197	3.4489	175.2224	0.4183	0.0583	G		
16	umc128.2	1	227602208	1	4.00E-04	4.54E-05	0.008	175.5028	0.4221	1.09E-08	A		
17	PZA03001.18	1	231037244	1	5.30E-05	4.54E-05	0.008	175.5028	0	2.35E-11	A		
18	PZA03001.15	1	231037305	1	8.00E-04	4.54E-05	0.008	175.5027	0.2966	2.82E-08	T		

Figure 3.2.12.Results in FASTmrEMMA module (**Result1**)

The **Result2** table with ten columns shows the final results of the FASTmrEMMA (a fast multi-locus random-SNP-effect EMMA) method. The corresponding column names are as follows: **Marker** (marker name or reference sequence number), **Chromosome**, **Marker position (bp)** in the chromosome, **QTN effect**, **LOD score**, **MAF** (minor allele frequency), **r² (%)** (the proportion of phenotypic variance explained by the putative QTN), **genotype for code 1**, **Var_Error** (residual error variance), and **Var_phen (total)** (total phenotypic variance), respectively.

	A	B	C	D	E	F	G	H	I	J	K
1	Marker	Chromosome	Marker position (bp)	QTN effect	LOD score	MAF	r ² (%)	Genotype for code 1	Var_Error	Var_phen (total)	
2	PZA03188.4	1	280719882	12.847	5.4526	0.2091	6.5419	C	343.0437	417.265	
3	PZA03559.1	2	15810363	8.3853	3.9294	0.3992	4.0025	C			
4											
5											
6											
7											
8											
9											
10											

Figure 3.2.13.Results in FASTmrEMMA module (**Result2**)

3.2.4 Manhattan plot

Click **Manhattan plot** button to preview Manhattan plot in an independent dialog window. Before saving the Figure, please set **the width and height of the Figure**, with the unit of pixel (px). And set **word resolution in the Figure**, with the unit of 1/72 inch, being pixels per inch (ppi). And set **figure resolution in the Figure**, with the unit of pixels per inch (ppi). The colors of the two adjacent chromosomes can be changed via the combo box, with a drop-down option. Set the critical value for $-\log_{10}(P)$, which is defaulted the value of $-\log(0.05/m_e)$, where m_e is the effective number of markers. Of course, users may change this value. Use **Save** button to choose a path and to save the

Figure, with three frequently used image formats: *.png, *.tiff and *.jpeg.




Figure 3.2.14. The Manhattan Plot of FASTmrEMMA module

3.2.5 QQ plot

Use **QQ plot** button to preview the QQ plot in an independent dialog window. The settings of the width, height, word resolution and figure resolution are the same as those in the Manhattan plot.




Figure 3.2.15. The QQ Plot of FASTmrEMMA module

Figure 3.2.14 and **Figure 3.2.15** show the preview picture of Manhattan plot and QQ plot with the parameter settings **in general resolution**, respectively.

If you want to obtain **high resolution** figure, we recommend the parameter settings shown in **Figure 3.2.16** and **Figure 3.2.17** for Manhattan plot and QQ plot, respectively. It is emphasized that preview figure may not seems same with **saved figure**. To give priority to saved figure in application.

Width (px):	<input type="text" value="10000"/>	Height (px):	<input type="text" value="6000"/>
Word resolution (1/72 inch, ppi):	<input type="text" value="30"/>	Figure resolution (ppi):	<input type="text" value="300"/>
Chromosome color (odd):	<input type="text" value="blue"/>	Chromosome color (even):	<input type="text" value="red"/>
Critical value for Manhattan Plot:	<input type="text" value="3"/>	<input type="button" value="Save"/> <input type="button" value="Cancel"/>	

Note: The parameter settings of preview figure is for general resolution, see manual for parameter settings of high resolution figure (*.png,*.tiff,*.jpeg).

Figure 3.2.16. Parameter settings for high resolution Manhattan plot of FASTmrEMMA Module

Width (px):	<input type="text" value="10000"/>	Height (px):	<input type="text" value="6000"/>
Word resolution (1/72 inch, ppi):	<input type="text" value="30"/>	Figure resolution (ppi):	<input type="text" value="300"/>
Point color:	<input type="text" value="black"/>	Line color:	<input type="text" value="red"/>
		<input type="button" value="Save"/> <input type="button" value="Cancel"/>	

Note: The parameter settings of preview figure is for general resolution, see manual for parameter settings of high resolution figure (*.png,*.tiff,*.jpeg).

Figure 3.2.17. Parameter settings for high resolution QQ plot of FASTmrEMMA module

3.3 GCIM module

Click the button **“Genome-wide Composite Interval Mapping (GCIM)”**, then the following dialog will appear.

Genome-wide Composite Interval Mapping (GCIM)

File Data_Format Graphic Help

Critical LOD score: Walk Speed for Genome-wide Scanning (cM):

Trait Selection:

About the software ⓘ

1. GCIM is a R software package for multi-locus linkage analysis based on Genome-wide Composite Interval Mapping.
2. Please cite: Wang Shi-Bo, Wen Yang-Jun, Ren Wen-Long, Ni Yuan-Li, Zhang Jin, Feng Jian-Ying, Zhang Yuan-Ming*. Mapping small-effect and linked quantitative trait loci for complex traits in backcross or DH populations via a multi-locus GWAS methodology. Scientific Reports 2016, 6: 29951.
3. The software package is developed by Yuan-Li Ni, Wen-Long Ren, Shi-Bo Wang & Yuan-Ming Zhang.

Version 2.0, Released Dec 2016

Figure 3.3.1. Screenshot of GCIM module GUI

3.3.1 Input data and Parameter setting

3.3.1.1 GCIM Format

To click GCIM will appear GCIM dialog box. In the dialog box, there are four steps. Firstly, the **Genotype**, **Phenotype** and **Linkage Map** buttons are used to input genotype, phenotype and linkage map datasets, respectively. Once one file is successfully uploaded, one tabbed page is added to the software notebook. Secondly, the user selects the required the population type, including **BC₁ ($F_1 \times P_1$)**, **BC₂ ($F_1 \times P_2$)**, **DH**, **RIL** and **Chromosome Segment Substitution Lines (CSSL)**. Thirdly, users should give the symbols in the marker translation table that users indicate genotypes. Users can enter any alpha numeric character(s) as these symbols. In this software, these symbols will translate into numeric variables in our method. Fourthly, the user selects the required **QTL-effect model type**: **Random** and **Fixed** models. Finally, some things will be implemented in this step: 1) to upload the above three files; 2) to select population type; 3) to conduct marker-genotype translation; and 4) to select QTL-effect model. Once users press the **DO** button, these things may be executed, until GCIM Format dialog box disappear and **info** window is appeared. When click **Continue** button, the next work is to select covariate and to run the program. If click **Cancel** button, the next work is not to execute these above things and GCIM Format dialog box is also disappeared.




Figure 3.3.2. GCIM Format Dataset dialog of GCIM module

3.3.1.2 QTL Ici Mapping Format

To click QTL Ici Mapping Format will appear QTL Ici Mapping dialog box. In the dialog box, there are three steps. Firstly, **QTL Ici Mapping Format** button is used to input QTL Ici Mapping format datasets. Once the file is successfully uploaded, one tabbed page is added to the software notebook. Secondly, the user selects the required population type (see §3.1.1). Thirdly, the users select the required QTL-effect model type (see §3.1.1). Finally, two things will be implemented in this step: 1) to upload QTL Ici Mapping-formatted files; and 2) to select population and model types. Once users press the **DO** button, the two things may be executed. Three tabbed pages (Genotype, Phenotype and Linkage map) will be added to the software notebook, and QTL Ici Mapping Format dialog box will be disappeared. To click **Cancel** button is not to execute the above things and QTL Ici Mapping Format dialog box is also disappeared.




Figure 3.3.3. QTL Ici Mapping Format Dataset dialog of GCIM module

3.3.1.3 WinQTL Cart Format

To click WinQTL Cart will appear WinQTL Cart dialog box. In the dialog box, there are three steps. Firstly, **WinQTL Cart Format** button is used to input WinQTL Cart format datasets. Once the related files are successfully uploaded, the tabbed pages are added to the software notebook. Secondly, the users select the required population type (see §3.1.1), and the users select the required model type (see §3.1.1). Finally, two things will be implemented in this step: 1) to upload WinQTL Cart format files;

and 2) to select population and model types. Once users press the **DO** button, the two things may be executed, three tabbed pages (Genotype, Phenotype and Linkage map) will be added to the software notebook, and WinQTLCart Format dialog box will be disappeared. If to click **Cancel** button, the next work is not to execute the above two things and WinQTLCart Format dialog box is also disappeared.




Figure 3.3.4. WinQTLCart Format Dataset dialog of GCIM module

3.3.1.4 Covariate

To click **Covariate** button will appear **Covariate** dialog box. In the dialog box, covariate(s) may be not included in the genetic model if there is no covariate. If not, the covariates should be included in the genetic model. Once the covariate file is successfully uploaded, the tabbed page will add to the software notebook.




Figure 3.3.5. The Covariate Format Dataset dialog of GCIM module

3.3.2 Run Program

Before run the program, please set up the **Critical LOD score**, which is defaulted by the 2.5 (LOD) value. If a LOD score for a putative QTL is larger than the Critical

LOD score, the putative QTL is viewed as true. If not, the putative QTL is viewed as false. Of course, the critical LOD score may be determined by permutation tests. If doing so, the software **WinQTLCart** is available. **Walk Speed for Genome-wide Scanning (cM)** is defaulted by the 1.0 (cM) value, which may be modified by users. **Trait Selection** is defaulted by the first trait. If user wants to analyze the third trait, the value in Trait Select column is set up by “3”. In the GCIM, to click Run button will execute the software. If the program runs, a progress bar with the “**Please be patient...**” words will appear in the bottom of the interface. If the program finished, a bar with the “**All done.**” will appear.




Figure 3.3.6. A running program interface of GCIM module




Figure 3.3.7. A finished program interface (the GCIM Results: [Trait 1](#))

3.3.3 Output results

To click **save** button in the **File** menu is used to save the result as *.csv file. When a trait is analyzed, the results may be saved, a dialog is used to select the pathway and the file name for the saved file.




Figure 3.3.8. The Save dialog of GCIM module

The **Results** are listed in ten columns for the GCIM (Genome-wide Composite Interval Mapping) method. The corresponding column names are as follows:

TraitID: Trait ID represented by an integer number.

Chr: Raw name of chromosome or chromosome ID represented by an integer number.

Position (cM): The scanning position (cM) on the chromosome.

Additive Effect: Estimated additive effect of the putative QTL.

LOD: LOD score for the putative QTL.

Left_Marker: Left flanking marker name of current scanning position (or putative QTL).

Right_Marker: Right flanking marker name of current scanning position (or putative QTL).

Var_Genet_(i): Genetic variance for all the detected QTL

r2 (%): proportion of phenotypic variance explained by single QTL.

Var_Error: residual variance for the full model.

Var_Phen (total): Phenotypic variance.

	A	B	C	D	E	F	G	H	I	J	K
1	TraitID	Chr	Position (cM)	Additive Effect	LOD	Left_Marker	Right_Marker	Var_Genet_(i)	r2 (%)	Var_Error	Var_Phen (total)
2	1	6	91.2693	5.4657	3.5896	*stm544tgagl	wPt-3728	29.8739	8.7301	123.7812	342.1955
3	1	11	49.4869	-12.1099	22.9343	cdo669b	csME1	146.6502	42.8557		
4	1	13	140.0322	-5.5446	6.1393	psr426	P40 N54-4	30.7425	8.9839		
5	1	5	64.1152	-3.3388	2.6002	mvg950	psr596	11.1477	3.2577		
6											
7											
8											
9											
10											
11											
12											
13											

Figure 3.3.9. The results of in the GCIM

3.3.4 Draw plot

If the running is ended, user may visualize the result by selecting **plot** button in the Graphic menu. Before clicking **Save** button in the **Genome-wide composite interval mapping (GCIM) figure** window, you can set parameter according to your demand. The parameter settings in Figure 3.3.10 are for **general resolution**.




Figure 3.3.10. The Genome-wide composite interval mapping (GCIM) figure using example data

Using example dataset in this software, the preview GCIM plots with the general resolution parameter settings were shown in **Figure 3.3.10**. Using real dataset of the Triticale DH population in Wurschum et al. (2014) Genetics, the preview GCIM plot with the general resolution parameter settings was shown in **Figure 3.3.11**.




Figure 3.3.11. The Genome-wide composite interval mapping (GCIM) figure using real data

If you want to obtain **high resolution** figure, we recommend the parameter settings shown in Figure 3.3.12. It is emphasized that preview figure may not seems same with **saved figure**. To give priority to save figure in application.

Legend and tick marks:	<input type="text" value="2.0"/>	LOD line size:	<input type="text" value="3.0"/>	Size for -log10(P) curve:	<input type="text" value="2.0"/>
Margin space:	<input type="text" value="2.5"/>	Space between tick marks and axis:	<input type="text" value="1.5"/>	Times for max(-log10(P)):	<input type="text" value="1.5"/>
Width (px):	<input type="text" value="10000"/>	Height (px):	<input type="text" value="4000"/>	Critical LOD score:	<input type="text" value="2.5"/>
Word resolution (1/72 inch, ppi):	<input type="text" value="12"/>	Figure resolution (ppi):	<input type="text" value="300"/>		
LOD line color:	<input type="text" value="red"/>	-log10(P) curve color:	<input type="text" value="gray50"/>	<input type="button" value="Save"/>	

Note: The parameter settings of preview figure is for general resolution, see manual for parameter settings of high resolution figure (*.png,*.tiff,*.jpeg).

Figure 3.3.12. Parameter settings for high resolution figure of GCIM module

4 REFERENCE

1. Wang Shi-Bo, Feng Jian-Ying, Ren Wen-Long, Huang Bo, Zhou Ling, Wen Yang-Jun, Zhang Jin, Jim M. Dunwell, Xu Shizhong*, Zhang Yuan-Ming*. Improving power and accuracy of genome-wide association studies via a multi-locus mixed linear model methodology. *Scientific Reports* 2016, **6**: 19444.
2. Wang Shi-Bo, Wen Yang-Jun, Ren Wen-Long, Ni Yuan-Li, Zhang Jin, Feng Jian-Ying, Zhang Yuan-Ming*. Mapping small-effect and linked quantitative trait loci for complex traits in backcross or DH populations via a multi-locus GWAS methodology. *Scientific Reports* 2016, **6**: 29951.
3. Wen Yang-Jun, Zhang Hanwen, Ni Yuan-Li, Huang Bo, Zhang Jin, Feng Jian-Ying, Wang Shi-Bo, Jim M. Dunwell, Zhang Yuan-Ming*, Wu Rongling*. Methodological implementation of mixed linear models in multi-locus genome-wide association studies. *Briefings in Bioinformatics* 2016, DOI: 10.1093/bib/bbw145